



TECHNICAL PROTOCOL

FOR

PGK-gb2-neo

**Pro- and Eukaryotic
Neomycin Selection Cassette**

(A001)

Gene Bridges GmbH
Im Neuenheimer Feld 584
69120 Heidelberg, Germany
Tel + 49 (0)6221 13708 11
Fax + 49 (0)6221 13708 29
Email: contact@genebridges.com
www.genebridges.com

CONTENTS

1 Eppendorf tubes + manual

1. PGK- gb2- neo: PCR template (50 ng/μl, 20μl)
2. This manual

Store tube at -20°C

Please read

The products listed in this manual are for research purposes only. They are not designed for diagnostic or therapeutic use in humans, animals or plants. The Red[®]/ET[®] recombination technology is the intellectual property of Gene Bridges GmbH.

Conditions of use

3.1 Purchaser will not manufacture, copy, reproduce, transmit, distribute, sell, lease, transfer, or improve upon the MATERIALS without prior written consent from GENE BRIDGES.

3.2 All MATERIALS relating Technologies shall be purchased from GENE BRIDGES or its authorized distributors. Use of any of the stated products from a source other than GENE BRIDGES will exempt GENE BRIDGES from any and all liabilities and warranties.

3.3 All MATERIALS purchased by research organizations, universities and other non-profit organizations may not be used for any commercial purpose. These MATERIALS are to be used for research purposes only. The MATERIALS may not be used to provide a commercial or non-commercial service, of any kind.

3.4 A purchase of MATERIALS by a private consumer is neither intended nor permitted.

Short Description:

“PGK-gb2-neo” template is designed to allow neomycin/kanamycin selection in prokaryotic and eukaryotic cells.

The PGK-gb2-neo template encodes the neomycin/kanamycin resistance gene which combines a prokaryotic promoter (gb2) for expression of kanamycin resistance in *E.coli* with a eukaryotic promoter (PGK) for expression of neomycin resistance in mammalian cells.

The prokaryotic promoter gb2 is a slightly modified version of the Em7 promoter; it mediates higher transcription efficiency than the normally used Tn5 promoter. The promoter of the mouse Phosphoglucokinase gene (PGK) is used as eukaryotic promoter. A synthetic polyadenylation signal terminates the kanamycin/neomycin expression.

Using the provided PCR template one can easily create a PGK-gb2-neo cassette flanked by any restriction sites to clone the cassette into the vector of choice. The restriction sites can be introduced by adding the corresponding sequence in the PCR primer. The template can easily be used to generate targeting constructs mediated by a single Red/ET Recombination step.

The “PGK-gb2-neo template” is not linear but plasmid based (3377bp in size). Due to its R6K origin it cannot replicate in most of the *E. coli* strains. The PCR product can therefore be used directly for downstream applications without any further purification.

At least 20 PCR reactions can be performed using 1µl per reaction as template.

Map: PGK-gb2-neo-template



NotI

```

1  GCGGCCGC ATTCTACCGG GTAGGGGAGG CGCTTTTCCC AAGGCAGTCT GGAGCATGCG CTTTAGCAGC CCCGCTGGGC
79  ACTTGCGCCT ACACAAGTGG CCTCTGGCTC GCACACATTG CACATCCACC GGTAGGCGCC AACCGGCTCC GTTCTTTGGT
159 GGGCCCTTCG CGCCACCTTC CACTCCTCCG CTAGTCAGGA AGTTCCTCCC CGCCCGCAG CTGCGTCTGT GCAGGACGTG
239 ACAAAATGAA GTAGCAGTTC TCACTAGTCT CGTGCAGATG GACAGCACCG CTGAGCAATG GAAGCGGGTA GGCCCTTTGGG
319 GCAGCGGCCA ATAGCAGCTT TGCTCCTTCG CTTTCTGGGC TCAGAGGCTG GGAAGGGGTG GGTCCGGGGG CGGGCTCAGG
399 GGC GGCTCA GGGCGGGGC GGGCGCCCGA AGGTCTCCG GAGGCCCGGC ATTCTGCACG CTTCAAAGC GCACGTCTGC
479 CGCGCTGTTT TCCTCTTCTT CATCTCCGGG CCTTTCGACC TGCAGC AGCAGGTGTT GACAATTAAT CATCGGCATA
555 GTATATCGGC ATAGTATAAT ACGACAAGGT GAGGAATAAA ACC ATG GGA TCG GCC ATT GAA CAA GAT GGA TTG
      i Met Gly Ser Ala Ile Glu Gln Asp Gly Leu
628 CAC GCA GGT TCT CCG GCC GCT TGG GTG GAG AGG CTA TTC GGC TAT GAC TGG GCA CAA CAG ACG ATC
111 His Ala Gly Ser Pro Ala Ala Trp Val Glu Arg Leu Phe Gly Tyr Asp Trp Ala Gln Gln Thr Ile
694 GGC TGC TCT GAT GCC GCC GTG TTC CCG CTG TCA CCG CAG GGG CGC CCG GTT CTT TTT GTC AAG ACC
333 Gly Oys Ser Asp Ala Ala Val Phe Arg Leu Ser Ala Gln Gly Arg Pro Val Leu Phe Val Lys Thr
760 GAC CTG TCC GGT GCC CTG AAT GAA CTG CAG GAC GAG GCA GCG CCG CTA TCG TGG CTG GCC ACG ACG
555 Asp Leu Ser Gly Ala Leu Asn Glu Leu Gln Asp Glu Ala Ala Arg Leu Ser Trp Leu Ala Thr Thr
826 GGC GIT CCT TGC GCA GCT GTG CTC GAC GTT GTC ACT GAA GCG GGA AGG GAC TGG CTG CTA ITG GGC
777 Gly Val Pro Oys Ala Ala Val Leu Asp Val Val Thr Glu Ala Gly Arg Asp Trp Leu Leu Leu Gly
892 GAA GTG CCG GGG CAG GAT CTC CTG TCA TCT CAC CTT GCT CCT GCC GAG AAA GTA TCC ATC ATG GCT
999 Glu Val Pro Gly Gln Asp Leu Leu Ser Ser His Leu Ala Pro Ala Glu Lys Val Ser Ile Met Ala
958 GAT GCA ATG CCG CGG CTG CAT ACG CTT GAT CCG GCT ACC TGC CCA TTC GAC CAC CAA GCG AAA CAT
1212 Asp Ala Met Arg Arg Leu His Thr Leu Asp Pro Ala Thr Oys Pro Phe Asp His Gln Ala Lys His
1024 CGC ATC GAG CGA GCA CGT ACT CCG ATG GAA GCC GGT CTT GTC GAT CAG GAT GAT CTG GAC GAA GAG
1433 Arg Ile Glu Arg Ala Arg Thr Arg Met Glu Ala Gly Leu Val Asp Gln Asp Asp Leu Asp Glu Glu
1090 CAT CAG GGG CTC CCG CCA GCC GAA CTG TTC GCC AGG CTC AAG GCG CGC ATG CCC GAC GGC GAG GAT
1655 His Gln Gly Leu Ala Pro Ala Glu Leu Phe Ala Arg Leu Lys Ala Arg Met Pro Asp Gly Glu Asp
1156 CTC GTC GTG ACC CAT GGC GAT GCC TGC TTG CCG AAT ATC ATG GTG GAA AAT GGC CGC TIT ICT GGA
1877 Leu Val Val Thr His Gly Asp Ala Cys Leu Pro Asn Ile Met Val Glu Asn Gly Arg Phe Ser Gly
1222 TTC ATC GAC TGT GGC CGG CTG GGT GTG CCG GAC CGC TAT CAG GAC ATA GCG TTG GCT ACC CGT GAT
2055 Phe Ile Asp Oys Gly Arg Leu Gly Val Ala Asp Arg Tyr Gln Asp Ile Ala Leu Ala Thr Arg Asp
1288 ATT GCT GAA GAG CTT GGC GGC GAA TGG CCT GAC CGC TTC CTC GTG CTT TAC GGT ATC GCC GCT CCC
2311 Ile Ala Glu Glu Leu Gly Gly Glu Trp Ala Asp Arg Phe Leu Val Leu Tyr Gly Ile Ala Ala Pro
1354 GAT TCG CAG CGC ATC GCC TTC TAT CGC CTT CTT GAC GAG TTC TTC TGA GCGGACTCTGGGGTTCGAATAAA
2533 Asp Ser Gln Arg Ile Ala Phe Tyr Arg Leu Leu Asp Glu Phe Phe
1426 GACCGACCAAGCGAC GTC TGA GAGCTCCCTG GCGAATTCGG TACCAATAAA AGAGCTTTAT TTTGATGATC
      XhoI
1497 TGTGTGTTGG TTTTGTGTG CGGCGCG CTCGAG

```

Please take into consideration that the sequence above does not reflect the complete plasmid but refers to the functional cassette.