



TECHNICAL PROTOCOL

FOR

**705-Cre**

expression plasmid

(A112)

Gene Bridges GmbH  
Im Neuenheimer Feld 584  
69120 Heidelberg, Germany  
Tel + 49 (0)6221 13708 11  
Fax + 49 (0)6221 13708 29  
Email: [contact@genebridges.com](mailto:contact@genebridges.com)  
[www.genebridges.com](http://www.genebridges.com)

## **CONTENTS**

### **1 Eppendorf tubes + manual**

1. 705-Cre: expression plasmid for Cre recombinase (0.2 µg/µl, 20 µl)
2. This manual

### **Store tube at -20°C**

### **Please read**

The products listed in this manual are for research purposes only. They are not designed for diagnostic or therapeutic use in humans, animals or plants. The Red<sup>®</sup>/ET<sup>®</sup> recombination technology is the intellectual property of Gene Bridges GmbH.

### **Conditions of use**

3.1 Purchaser will not manufacture, copy, reproduce, transmit, distribute, sell, lease, transfer, or improve upon the MATERIALS without prior written consent from GENE BRIDGES.

3.2 All MATERIALS relating Technologies shall be purchased from GENE BRIDGES or its authorized distributors. Use of any of the stated products from a source other than GENE BRIDGES will exempt GENE BRIDGES from any and all liabilities and warranties.

3.3 All MATERIALS purchased by research organizations, universities and other non-profit organizations may not be used for any commercial purpose. These MATERIALS are to be used for research purposes only. The MATERIALS may not be used to provide a commercial or non-commercial service, of any kind.

3.4 A purchase of MATERIALS by a private consumer is neither intended nor permitted.

### **Short Description:**

705-Cre plasmid is designed for use in Cre-mediated genomic manipulations. The plasmid has a pSC101 origin which maintains low copy and replicates at 30°C. The plasmids will not propagate and will get lost when incubated at 37°C.

The expression of the Cre-recombinase is driven by the thermosensitive promoter cl578 ( $\lambda_{PR}$  promoter). Therefore, the expression of Cre is repressed at 30°C and induced between 37-42°C.

The plasmid carries a chloramphenicol resistance.

### **Note:**

The sequence of 705-Cre was compiled from information found in the sequence databases, published literature, and other sources, together with partial sequences obtained by Genebridges. This vector has not been completely sequenced.

The digestion patterns for *Bam*HI, *Eco*RI, *Hind*III, *Pst*I, *Xba*I and *Xho*I are indicated

### **Reference:**

Buchholz, F., Angrand, P.-O. and Stewart, A.F. (1996) „A simple assay to determine the functionality of Cre or FLP recombination targets in genomic manipulation constructs” *Nucleic Acids Research* 24, 3118-3119.


Zhang, Y., Buchholz, F., Muyrers, J.P.P. and Stewart, A.F. (1998) “A new logic for DNA engineering using recombination in *Escherichia coli*” *Nature Genetics* 20, 123-128.

## Site Specific Recombination to Remove Selection Marker.

1. 705-Cre plasmid is transformed into an *E.coli* strain, which contains a targeting plasmid carrying a floxed selection marker (e.g. pPGK-neo resistance gene).
2. After transformation (electroporation or heat shock), add 1 ml of LB medium to the tube and incubate at 30°C for 1.5 hr with shaking.
3. Streak out the cells on L.B. plates containing 15 µg/ml of chloramphenicol (cm) plus ampicillin (amp; selection marker for the targeting plasmid).
4. Incubate at 30°C for more than 24 hours (since the colonies grow slowly).
5. Pick a single colony and grow the cells in 1 ml of LB medium with 50 µg/ml of amp (resistance of the targeting plasmid) at 30°C for 2-3 hours.
6. Switch temperature to 37°C and incubate overnight.  
  
(During incubation at 37°C, Cre protein is expressed and the loxP sites recombined, at the same time, 705-Cre plasmid is lost.)
7. Prepare plasmid DNA and digest part of the DNA to check the restriction pattern
8. Re-transform the checked DNA to remove the unrecombined plasmid.

About 95% of floxed fragment will be recombined. Step 8 is therefore important to obtain the pure and recombined plasmid.


**Map:**




### Restriction pattern plasmid 705-Cre, 705-Flp, 706-Cre and 706-Flp

[bp]	705-Cre	705-Flp	706-Cre	706-Flp
<i>Bam</i> HI	9700	undigested	1300 4200 5800	11000
<i>Eco</i> RI	1100 4300 (2x)	1800 4300 (2x)	1000 4200 6100	1700 4000 5300
<i>Hind</i> III	1100 1600 7000	1300 1800 7100	1000 1300 1600 3200 4200	1200 1600 3600 4600
<i>Pst</i> I	1300 8400	10200	1300 10000	11000
<i>Xba</i> I	undigested	undigested	undigested	undigested
<i>Xho</i> I	undigested	undigested	undigested	undigested


#### 705-Cre




#### 705-Flp




#### 706-Cre



#### 706-Flp



### Functional test 705-Cre, 705-Flp, 706-Cre and 706-Flp



The functional test was performed as described on page 4 of the manual. The plasmids pSVpaZ11 (size 7.3 kb; with a 1.1 kb FRT flanked fragment) and pSVpaX1 (size 7.3kb; with a floxed 1.1 kb fragment) were used as targeting plasmids. Miniprep. DNA from two colonies was isolated and the targeting plasmids linearized by *NotI* digestion to check for successful recombination.

**705-Cre, 706-Cre:** The size of pSVpaX1 shows the successful recombined size of 6.2 kb (lanes 1+2, lanes 3+4, respectively). The negative control pSVpaX1 shows the original 7.3kb band (lane 9).

**705-Flp, 706-Flp:** The size of pSVpaZ11 shows the successful recombined size of 6.2 kb (lane 5+6, lanes 7+8, respectively). pSVpaZ11 is used as negative control and shows the original size of 7.3 kb (lane 10).